

Baňka Erlenmeyerova se zábrusem | GLASSCO

Kónická baňka se zábrusem bez zátky

- Bílá stupnice a popisovací plocha
- Dle normy ISO 4797

Technické parametry

Materiál: borosilikátové sklo 3.3



Baňka Erlenmeyerova se zábrusem

Objem (ml)	Zábrus NZ	Průměr (mm)	Výška (mm)	Balení (ks)	Katalogové číslo
25	14/23	42	70	1	N120206
25	19/26	42	70	1	N120207
50	14/23	51	85	1	N120208
50	19/26	51	85	1	N120209
50	29/32	51	85	1	N120210
100	14/23	64	105	1	N120211
100	19/26	64	105	1	N120212
100	24/29	64	105	1	N120213
200	29/32	79	130	1	N120217
250	29/32	85	140	1	N120220
300	29/32	87	155	1	N120221
500	29/32	105	175	1	N120223
1 000	29/32	131	215	1	N120226
2 000	29/32	166	275	1	N120228
3000	45/40	188	310	1	N120229
5000	45/40	220	365	1	N120230